

## Features at a glance:

- 2 years warranty
- durable, rugged design
- simple operation
- low noise

# oxy 5000

## Oxygen Concentrator

Made in Germany



*silent operation*  
*< 40 db(A)*

*modular*  
*component*  
*setup*

*USB interface*

*microprocessor controlled*

*easy-to-operate*

*with oxygen sensor*



### Shipping contents

oxy 5000 O<sub>2</sub> concentrator

O/N 48-1000

including:

concentrator 48-1001

humidifier 48-1502

nasal cannula 2m 41-1506

nasal cannula 5m 41-1507

instructions for use 48-1201

*more accessories on request*

### Specifications

dimensions: 600 x 290 x 400 mm

weight: 15 kg

power: 230 V / 50 Hz

power consumption: 295 W

fuses, primary: T 3.15 AH 250 V

fuses, secondary: T 1 AL 250 V

operating temperature: +10° - +40° C

storage temperature: -25° - +70° C

sound power: < 40 dB(A)

MDD Classification: IIa

flow: 0 to 5 lpm

oxygen: 1 to 4 lpm: 95 % - 3 %

3 to 4 lpm: 93 % +/- 3 %

4 to 5 lpm: 85 % +/- 3 %

alarms: temperature, power failure

oxygen, sensor, system

dust filter: in the back

fine filter: in the back

hours counter: LCD

interface: USB

warranty: 2 years

Bitmos GmbH  
Himmelgeister Str. 37  
D 40225 Düsseldorf Germany

Phone: +49 211 60101030  
Fax: +49 211 60101050

[sales@bitmos.de](mailto:sales@bitmos.de)  
[www.bitmos.de](http://www.bitmos.de)

**Bitmos**  
MEDIZINTECHNIK